ENGLISH WARNINGS

The Muscle Massager is intended for commercial or home use. We strive to make the Muscle Massager as safe as possible. This is an advanced mechanical device with electric components. If the Muscle Massager and its accessories are not used or maintained properly, there is a risk of fire, electric shock, or injury. Do not place your finger or any objects near the metal piston above the attachment while the Muscle Massager is in use or immediately after use. The metal piston can become extremely hot. Do not use any attachments other than the ones provided by Flexir Recovery and only as instructed by Flexir Recovery. Do not use the Muscle Massager on your head or near your genitals. Do not use the Muscle Massager above the waist during pregnancy. Do not get hair caught in any moving part of the Muscle Massager. Do not use the Muscle Massager above your Adam's apple or C4. Do not get the Muscle Massager wet. Only clean the Muscle Massager with a magic eraser or a lightly damp towel. Do not block the vents of the motor. Never attempt to disassemble the Muscle Massager. Never store the Muscle Massager at temperatures above 40°C/104°F. Only charge the Muscle Massager when the ambient temperature is between 0°C/32°F and 40°C/104°F. Never expose the Muscle Massager to water. For long-term storage, the Muscle Massager should be fully charged. This product is not intended to diagnose, treat, or prevent any diseases. BEFORE USING OR CHARGING THE Muscle Massager, READ ALL THE INSTRUCTIONS AND CAUTIONARY MARKINGS IN THIS MANUAL AND ON THE Muscle Massager. When using the Muscle Massager, these basic precautions should always be followed including:

1. USE ONLY AS INSTRUCTED. Use the Muscle Massager as described in this Muscle Massager Instruction Manual. Do not carry out any maintenance other than what is shown in the Muscle Massager Instruction Manual or advised by Flexir Recovery.
2. NOT FOR CHILDREN. The Muscle Massager is not intended for use by young children or persons with reduced physical, sensory or reasoning capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction by a responsible person. Do not allow the Muscle Massager to be used as a toy. Close attention is necessary when used by or near children. Children should be supervised to ensure

that they do not play with the Muscle Massager.

3. DO NOT PULL OUT THE BATTERY WHILE THE MUSCLE MASSAGER IS ON.

4. REDUCING RISK OF INJURY. Charge the Muscle Massager by plugging cable included in the package on your USB AC adapter of choice. Muscle Massager can also be charged with a regular USB-C charger certified by a competent body. Do not use an AC adapter that is not certified or that has received a sharp blow, been dropped, or damaged in any way.

5. USE THE MUSCLE MASSAGER ONLY WITH FLEXIR RECOVERY MUSCLE MASSAGER BATTERIES. Use of other batteries may result in a risk of fire, electric shock, or personal injury.

6. CHARGING LOCATIONS. The Muscle Massager should be charged indoors in a well ventilated, dry location. Do not charge the Muscle Massager outdoors, in a bathroom or within 10 feet (3.1 meters) of a bathtub, or pool. Do not use the Muscle Massager on wet surfaces and do not expose the device to moisture, rain, or snow. Do not use the Muscle Massager in the presence of explosive atmospheres (gaseous fumes, dust, or flammable materials). Sparks may be generated, possibly causing a fire.

7. DO NOT OVERCHARGE. Do not leave the Muscle Massager connected to the adapter for more than an hour after the Muscle Massager has been fully charged. The Muscle Massager includes a protection system to avoid the risk of overcharging. However, overcharging the Muscle Massager may reduce its life over time.

8. DO NOT BURN OR INCINERATE THE Muscle Massager. The internal battery may explode, causing personal injury or damage. Toxic fumes and materials are created when the battery is burned.

9. DO NOT CRUSH, DROP, OR DAMAGE THE Muscle Massager charging cable. Do not use a charging cable that has received a sharp blow, been dropped, run over, or damaged in any way.

10. BATTERY CHEMICALS CAUSE SERIOUS BURNS. Never allow the internal battery to come into contact with the skin, eyes, or mouth. If a damaged battery leaks chemicals, use rubber or neoprene gloves to dispose of it. If skin is exposed to battery fluids, wash with soap and water then rinse with vinegar. If eyes are exposed to battery chemicals, immediately rinse with water for 20 minutes and seek medical attention. Remove and dispose of contaminated clothing.

11. DO NOT OPERATE UNDER BLANKET AND PILLOW. Excessive heating can occur and cause fire, electric shock, or injury to persons.
12. STORING THE Muscle Massager. Store in a cool, dry place. Do not store the Muscle Massager where temperatures may exceed 40°C/104°F such as in direct sunlight, in a vehicle, or metal building during the summer.
13. DISPOSING OF BATTERIES. The internal lithium-ion batteries for the Muscle Massager are more environmentally friendly than some other types of batteries (e.g, nickel-cadmium). Always dispose of the battery according to government, regional, and local regulations. Contact a recycling agency in your area for recycling locations. Even discharged batteries contain some energy. Before disposing, use electrical tape to cover the Muscle Massager's 12V connector to prevent the battery from shorting, which could cause a fire or explosion.
14. DO NOT DISASSEMBLE. Disassembly or incorrect reassembly may result in the risk of electric shock, fire, or exposure to battery chemicals. The warranty will be void if the Muscle Massager is disassembled or if any parts have been removed.
15. SERVICE. If the Muscle Massager is not working properly, has received a sharp blow, has been dropped, damaged, left outdoors, or dropped into water, do not use. Do not attempt to repair or disassemble the Muscle Massager, which may result in an electric shock or fire.
16. DO NOT USE WHILE BATHING OR IN THE SHOWER, TUB, OR SINK. Do not place or store the Muscle Massager where it can fall or be pulled into a tub or sink. Do not place in or drop into water or other liquid. Do not reach for an appliance that has fallen into water. Unplug immediately.
17. THERMAL LIMITER. The Muscle Massager has an automatic reset thermal limiter that shuts off the Muscle Massager to prevent overheating and fire.

This device can be used by children aged from 8 years and above, and persons with reduced physical, sensory or mental capabilities, or lack of experience and knowledge if they have been given supervision or instruction concerning use of the device in a safe way and understand the hazards involved. Children shall not play with the device. Cleaning and user maintenance shall not be used without supervision. This device contains batteries that are only replaceable by authorized personnel. This device contains batteries that are non-replaceable.

LIMITED WARRANTY

For full warranty information, please visit FlexirRecovery.com/warranty. To request a copy of the warranty by mail, you may send a request to the following address:

Flexir Recovery
Attn: Warranty Department
251 Little Falls Drive
Wilmington, DE 19807

Please note, this is not a return address or a retail location. No products or packages will be accepted at this location.

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Flexir, LLC. The Masonry, 151-156 Thomas St Suite 208, Dublin, Ireland